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### *What should be done?*

- *Do not deliver milk with a high quantity of cells.*
- *Do the milk analysed by the dairy laboratory and ask for the result.*
- *Veterinary treatment*
- *In case of infection, apart from looking for pathogens, causes could also be due to poor maintenance (cleaning, disinfection) of the milking machine.*

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*Allow regular checking of the milking machine by your dealer.*

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### Meaning of the quantities of cells in raw milk

*The milk cell rate determines the health condition of the udder. The constant regeneration and reaction in the gland require a certain (normal) quantity of cells in milk. However, it is not constant, but varies in the quarters between 10 000 and 200,000 cells per ml of milk.*

#### **Sick wards :**

- *Production of milk with high quantity of leukocyte cells.*
- *Less milk.*
- *Loss of income for the breeder.*

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### Why this test ?

*This test gives informations on the status of the quarters well before the change of aspect of the milk is visible to the naked eye (inflammation of the breast tissue, lump, water + blood, etc...). One can notice an early disease of the udders from 200,000 cells onwards. Regular monitoring will indicate a possible rapid increase in cells (several million).*

***It is therefore very important to detect the disease at the very beginning in order to prevent more damage.***

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### When should you use the test ?

- *Once every 4 weeks for prevention*
- *After the result from the laboratory dairy analysis in order to quickly detect sick animals*
- *Ten days after a treatment of mastitis to check the state of healing*
- *At the purchase of an animal*
- *Two weeks after calving*

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### How to use the test ? *The test should be done just before milking.*

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*Important: Remove the first jets, as the initial milk contains a high rate of cells.*

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*Pull from each quarter a little of milk in each bowl of the test tray identifying which quarter correspond to which bowl.*

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
*By inclination of the plate, the surplus of milk will be discharged to the scale.*

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*Prime properly the dosing bottle so that the necessary amount of test liquid is versed into each bowl.*

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*Throw the mixture and rinse the tray with clear water. Without wiping it, it is ready again for a second test.*

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*The mixture remains liquid (no lumps)*

*Evaluation : up to 200,000 cells : milk fulfils regulations.*

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*The mixture shows, by inclination of the tray, light coagulation that one clearly sees on the scale mark.*

*Evaluation : the quantity of cells increased slightly. 200,000 to 500,000 cells : beginning of probable infection.*

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*The mixture shows a high coagulation (gelatinous).*

*Evaluation : the quantity of cells increased sharply. 500,000 to 1,000,000 cells : mastitis.*

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*The mixture shows an exceptionally strong and gelatinous coagulation. In addition, can change colour to red/blue.*

*Evaluation : 1,000,000 to several million cells : very strong udder infection.*

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